

Use Case: Record Session

Primary Actor: Device User / Device Control

Scope: Oasis Pro System

Level: Subfunction

Stakeholders and interests:

- User of the device - to record their therapy session and have it saved to their user profile treatment history
- Oasis Pro Device - to gather the session type, duration, and intensity level and save it to the active users treatment history

Main success scenario:

- Control system receives signal from device button to record the therapy session
- User initializes the therapy session
- Control system carries out the therapy session
- Control system prompts the user if they wish to save their recorded session upon termination of the session
- Control system receives signal from user prompt that the user wishes to save their recorded session
- Control system saves the session information (session type, duration, and intensity level) to the active users treatment history
- Control system returns to the main session screen

Variations:

5'. Control system receives signal from user prompt that the user wishes to delete their recorded session - the control system will dismiss the recorded session and return to the main session screen

